



Certificate of Analysis

Sample: TE31019007-001

Harvest/Lot ID: 2308023B

Batch#: 2308023B

Batch Date: 10/19/23

Sample Size Received: 51.75 gram

Total Amount: 1 units

Retail Product Size: 300 gram

Ordered: 10/17/23

Sampled: 10/19/23

Completed: 10/25/23

Revision Date: 10/26/23

Oct 26, 2023 | e2e Pharma

3279 E. Harbour Drive
 Phoenix, AZ, 85034, US



PASSED

Pages 1 of 6


PRODUCT IMAGE



SAFETY RESULTS



Pesticides
PASSED



Heavy Metals
PASSED



Microbials
PASSED



Mycotoxins
PASSED



Residuals Solvents
PASSED



Filtration
PASSED



Water Activity
NOT TESTED



Moisture
NOT TESTED



Terpenes
NOT TESTED

MISC.

750 mg CBD



Cannabinoid

PASSED



Total THC
ND

Total THC/Container : 0.000 mg



Total CBD
0.2678%

Total CBD/Container : 803.400 mg



Total Cannabinoids
0.2932%

Total Cannabinoids/Container : 879.600 mg

	D9-THC	THCA	CBD	CBDA	CBG	CBGA	CBN	D8-THC	CBDV	THCV	CBC
%	ND	ND	0.2678	ND	0.0192	ND	0.0062	ND	<0.0070	ND	<0.0070
mg/unit	ND	ND	803.400	ND	57.600	ND	18.600	ND	<21.000	ND	<21.000
LOD	0.0020	0.0020	0.0020	0.0020	0.0020	0.0010	0.0010	0.0020	0.0020	0.0020	0.0010
	%	%	%	%	%	%	%	%	%	%	%

Analyzed by:
121, 272, 104

Weight:
1.0344g

Extraction date:
10/20/23 18:31:22

Extracted by:
121

Analysis Method : SOP.T.30.500, SOP.T.30.031, SOP.T.40.031

Analytical Batch : TE002937POT

Instrument Used : TE-245 "Muad'Dib" (Infused)

Analyzed Date : 10/20/23 18:57:08

Reviewed On : 10/25/23 10:53:04

Batch Date : 10/20/23 13:59:45

Dilution : 40

Reagent : 083123.21; 101723.R07; 101123.R01; 100623.R10; 072522.R32

Consumables : 947.100; 20220108; 8000031463; 00331867-5; 12541-224CD-224; 210630-306-D; 210725-598-D; 0000185478; GD220011

Pipette : TE-059 SN:20A04528 (20-200uL); TE-065 SN:20B18327 (100-1000uL); TE-164 SN: 21H24198 (Isopropanol)

Full spectrum cannabinoid analysis utilizing High Performance Liquid Chromatography with Photo Diode Array detector (HPLC-PDA) for analysis. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.031 for sample prep, SOP.T.40.031 for analysis on Shimadzu LC-20X0 series HPLCs). Potency results for cannabis flower products are reported on an "as received" basis, without moisture correction.

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors.

Ariel Gonzales
Lab Director

State License #
00000024LCMD66604568
ISO 17025 Accreditation # 97164



Signature
10/25/23

Revision: #1 - New photo



Certificate of Analysis

PASSED

e2e Pharma

Sample : TE31019007-001

Harvest/Lot ID: 2308023B

Batch# : 2308023B

Sampled : 10/19/23

Ordered : 10/19/23

Sample Size Received : 51.75 gram

Total Amount : 1 units

Completed : 10/25/23 Expires: 10/26/24

Sample Method : SOP Client Method

3279 E. Harbour Drive
Phoenix, AZ, 85034, US
Telephone: (602) 737-0077
Email: shannon.bard@e2epharmamfg.com

Page 2 of 6



Pesticides

PASSED

Pesticide	LOD	Units	Action Level	Pass/Fail	Result	Pesticide	LOD	Units	Action Level	Pass/Fail	Result
AVERMECTINS (ABAMECTIN B1A)	0.0170	ppm	0.5	PASS	ND	TOTAL SPINOSAD	0.0060	ppm	0.2	PASS	ND
ACEPHATE	0.0100	ppm	0.4	PASS	ND	SPIROMESIFEN	0.0080	ppm	0.2	PASS	ND
ACETAMIPRID	0.0050	ppm	0.2	PASS	ND	SPIROTETRAMAT	0.0060	ppm	0.2	PASS	ND
ALDICARB	0.0140	ppm	0.4	PASS	ND	SPIROXAMINE	0.0040	ppm	0.4	PASS	ND
AZOXYSTROBIN	0.0050	ppm	0.2	PASS	ND	TEBUCONAZOLE	0.0040	ppm	0.4	PASS	ND
BIFENAZATE	0.0060	ppm	0.2	PASS	ND	THIACLOPRID	0.0060	ppm	0.2	PASS	ND
BIFENTHRIN	0.0050	ppm	0.2	PASS	ND	THIAMETHOXAM	0.0060	ppm	0.2	PASS	ND
BOSCALID	0.0050	ppm	0.4	PASS	ND	TRIFLOXYSTROBIN	0.0060	ppm	0.2	PASS	ND
CARBARYL	0.0080	ppm	0.2	PASS	ND	CHLORFENAPYR *	0.0270	ppm	1	PASS	ND
CARBOFURAN	0.0050	ppm	0.2	PASS	ND	CYFLUTHRIN *	0.0150	ppm	1	PASS	ND
CHLORANTRANILIPROLE	0.0110	ppm	0.2	PASS	ND	Analyzed by: 152, 272, 104 Weight: 0.4957g Extraction date: 10/24/23 17:48:15 Extracted by: 152 Analysis Method: SOP.T.30.500, SOP.T.30.104.AZ, SOP.T.40.104.AZ Analytical Batch: TE002961PES Reviewed On: 10/25/23 15:34:10 Instrument Used: TE-118 *MS/MS Pest/Myco 1*, TE-261 *UHPLC - Pest/Myco 2* Batch Date: 10/23/23 15:47:49 Analyzed Date: 10/24/23 19:12:24 Dilution: 25 Reagent: 102023.R01; 091323.R20; 102323.R16; 101823.R02; 101123.R02; 063123.R02; 041823.06 Consumables: 947.100; 00334958-5; 00332484-2; 1008439554; 28521042; 210823-1124; 090623; 269336; M0040957; GD220011 Pipette: TE-056 SN:21D58687; TE-060 SN:20C35457 (20-200uL); TE-108 SN:20B18337 (100-1000uL) Pesticide screening is carried out using LC-MS/MS supplemented by GC-MS/MS for volatile pesticides. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.104.AZ for sample prep, and SOP.T.40.104.AZ for analysis on ThermoScientific Altis TSQ with Vanquish UHPLC).					
FLONICAMID	0.0090	ppm	1	PASS	ND	Analyzed by: 152, 272, 104 Weight: 0.4957g Extraction date: 10/24/23 17:48:15 Extracted by: 152 Analysis Method: SOP.T.30.500, SOP.T.30.104.AZ, SOP.T.40.154.AZ Analytical Batch: TE002975VOL Reviewed On: 10/25/23 15:34:41 Instrument Used: TE-091 *GC - Volatile Pesticides 1*, TE-094 *MS/MS - Volatile Pesticides 1* Batch Date: 10/24/23 17:50:08 Analyzed Date: 10/24/23 18:57:46 Dilution: 25 Reagent: 102023.R01; 091323.R20; 111921.03; 030623.03 Consumables: 947.100; 00334958-5; 00332484-2; 1008439554; 28521042; 210823-1124; 090623; 269336; M0040957; GD220011 Pipette: TE-056 SN:21D58687; TE-060 SN:20C35457 (20-200uL); TE-108 SN:20B18337 (100-1000uL) Supplemental pesticide screening using GC-MS/MS to quantitatively screen for Chlorfenapyr, Cyfluthrin, Cypermethrin, and Diazinon; as well as the qualitative confirmation of Dichlorvos, Permethrin, Piperonyl Butoxide, Prallethrin, Propiconazole, Pyrethrin, and Tebuconazole which are all quantitatively screened using LC-MS/MS. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.104.AZ for sample prep, and SOP.T.40.154.AZ for analysis using a ThermoScientific 1310-series GC equipped with a TriPlus RSH autosampler and detected on a TSQ 9000-series mass spectrometer).					
FLUDIOXONIL	0.0060	ppm	0.4	PASS	ND						
HEXYTHIAZOX	0.0050	ppm	1	PASS	ND						
IMAZALIL	0.0110	ppm	0.2	PASS	ND						
IMIDACLOPRID	0.0080	ppm	0.4	PASS	ND						
KRESOXIM-METHYL	0.0070	ppm	0.4	PASS	ND						
MALATHION	0.0070	ppm	0.2	PASS	ND						
METALAXYL	0.0040	ppm	0.2	PASS	ND						
METHIOCARB	0.0040	ppm	0.2	PASS	ND						
METHOMYL	0.0050	ppm	0.4	PASS	ND						
MYCLOBUTANIL	0.0100	ppm	0.2	PASS	ND						
NALED	0.0070	ppm	0.5	PASS	ND						
OXAMYL	0.0080	ppm	1	PASS	ND						
PACLOBUTRAZOL	0.0050	ppm	0.4	PASS	ND						
TOTAL PERMETHRINS	0.0030	ppm	0.2	PASS	ND						
PHOSMET	0.0100	ppm	0.2	PASS	ND						
PIPERONYL BUTOXIDE	0.0050	ppm	2	PASS	ND						
PRALLETHRIN	0.0130	ppm	0.2	PASS	ND						
PROPICONAZOLE	0.0050	ppm	0.4	PASS	ND						
PROPOXUR	0.0050	ppm	0.2	PASS	ND						
TOTAL PYRETHRINS	0.0010	ppm	1	PASS	ND						
PYRIDABEN	0.0040	ppm	0.2	PASS	ND						

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation, Limit of Detection (LOD) and Limit of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors.

Ariel Gonzales
Lab Director

State License #
00000024LCMD66604568
ISO 17025 Accreditation # 97164



Signature
10/25/23

Revision: #1 - New photo



Certificate of Analysis

PASSED

e2e Pharma

 3279 E. Harbour Drive
 Phoenix, AZ, 85034, US
 Telephone: (602) 737-0077
 Email: shannon.bard@e2epharmamfg.com

Sample : TE31019007-001

Harvest/Lot ID: 2308023B

Batch# : 2308023B

Sampled : 10/19/23

Ordered : 10/19/23


Sample Size Received : 51.75 gram

Total Amount : 1 units

Completed : 10/25/23 Expires: 10/26/24

Sample Method : SOP Client Method

Page 3 of 6



Residual Solvents

PASSED

Solvents	LOD	Units	Action Level	Pass/Fail	Result
BUTANES	168.2000	ppm	5000	PASS	ND
METHANOL	87.7000	ppm	3000	PASS	ND
PENTANES	163.9000	ppm	5000	PASS	ND
ETHANOL	142.2000	ppm		TESTED	ND
ETHYL ETHER	193.1000	ppm	5000	PASS	ND
ACETONE	37.6000	ppm	1000	PASS	ND
2-PROPANOL	156.2000	ppm	5000	PASS	ND
ACETONITRILE	12.2000	ppm	410	PASS	ND
DICHLOROMETHANE	22.7000	ppm	600	PASS	ND
HEXANES	8.4000	ppm	290	PASS	ND
ETHYL ACETATE	179.0000	ppm	5000	PASS	ND
CHLOROFORM	2.4100	ppm	60	PASS	ND
BENZENE	0.1150	ppm	2	PASS	ND
ISOPROPYL ACETATE	168.6000	ppm	5000	PASS	ND
HEPTANE	152.8000	ppm	5000	PASS	ND
TOLUENE	26.2000	ppm	890	PASS	ND
XYLENES	53.2000	ppm	2170	PASS	ND

Analyzed by: 93, 30, 104	Weight: 0.0247g	Extraction date: 10/20/23 18:51:05	Extracted by: 93
-----------------------------	--------------------	---------------------------------------	---------------------

Analysis Method : SOP.T.40.044.AZ

Analytical Batch : TE002934SOL

Instrument Used : TE-092 "GC - Solvents 1", TE-095 "MS - Solvents 1", TE-098 "Injector - Solvents 1", TE-100 "HS - Solvents 1", TE-113 "Vacuum Pump - Solvents 1"

Reviewed On : 10/23/23 12:03:09

Batch Date : 10/20/23 13:10:50

Analyzed Date : 10/20/23 15:49:58

Dilution : N/A

Reagent : 072722.01; 051223.03; 032023.03

Consumables : 428251; 19000-1; GD220011

Pipette : N/A

Residual solvents screening is performed using GC-MS which can detect below single digit ppm concentrations. (Method: SOP.T.40.044.AZ for sample prep and analysis via ThermoScientific 1310-series GC equipped with a TriPlus 500 Headspace autosampler and detection carried out by ISQ7000-series mass spectrometer). Butanes are reported as the sum of n-Butane and Isobutane. Pentanes are reported as the sum of n-Pentane, Isopentane, and Neopentane. Hexanes are reported as the sum of n-Hexane, 2-Methylpentane, 3-Methylpentane, 2,2-Dimethylbutane, and 2,3-Dimethylbutane. Xylenes are reported as the sum of Ethyl Benzene, m-Xylene, p-Xylene, and o-Xylene.



Certificate of Analysis

PASSED

e2e Pharma

3279 E. Harbour Drive
Phoenix, AZ, 85034, US
Telephone: (602) 737-0077
Email: shannon.bard@e2epharmamfg.com

Sample : TE31019007-001
Harvest/Lot ID: 2308023B

Batch# : 2308023B
Sampled : 10/19/23
Ordered : 10/19/23

Sample Size Received : 51.75 gram
Total Amount : 1 units
Completed : 10/25/23 Expires: 10/26/24
Sample Method : SOP Client Method

Page 4 of 6

	Microbial	PASSED		Mycotoxins	PASSED
---	------------------	---------------	---	-------------------	---------------

Analyte	LOD	Units	Result	Pass / Fail	Action Level
SALMONELLA SPP			Not Present in 1g	PASS	
ESCHERICHIA COLI REC	10.0000	CFU/g	ND	PASS	100
TYM	1.0000	Colonies	ND	TESTED	

Analyzed by: 96, 87, 104 **Weight:** 0.945g **Extraction date:** 10/20/23 14:05:11 **Extracted by:** 87,96
Analysis Method : SOP.T.40.056B, SOP.T.40.058.FL, SOP.T.40.208, SOP.T.40.209.AZ
Analytical Batch : TE002931MIC **Reviewed On :** 10/24/23 14:30:14
Instrument Used : TE-234 "bioMerieux GENE-UP" **Batch Date :** 10/20/23 12:21:49
Analyzed Date : 10/24/23 11:12:37
Dilution : 10
Reagent : 091123.07; 091123.10; 101323.R02; 101323.R03; 101923.04; 101923.05; 080423.33; 080423.37; 051623.41; 080423.42; 080323.03; 101823.04; 101823.05; 091323.20; 091323.21; 051923.17; 090723.04; 090723.05
Consumables : 22507; 33T5N9; 41310-229C4-229I; 1008439554; 211108-071-B; 210715-071; 28521042; 12541-224CD-224; 111521CH02; 210823-1124; 269336; X0028AKTV1; X002E5BZFT; 40172
Pipette : TE-053 SN:20E78952; TE-058 SN:20C35427; TE-061 SN:20C35454; TE-062 SN:20C50491; TE-066 SN:20D56970; TE-069 SN:21B23920; TE-109 SN:20B18330; TE-256 Dispensette S Bottle Top Dispenser SN:20G36073; TE-258

Analyte	LOD	Units	Result	Pass / Fail	Action Level
TOTAL AFLATOXINS	1.4870	ppb	ND	PASS	20
AFLATOXIN B1	1.4700	ppb	ND	PASS	20
AFLATOXIN B2	1.8000	ppb	ND	PASS	20
AFLATOXIN G1	1.9000	ppb	ND	PASS	20
AFLATOXIN G2	3.2500	ppb	ND	PASS	20
OCHRATOXIN A	4.6100	ppb	ND	PASS	20

Analyzed by: 152, 272, 104 **Weight:** 0.4957g **Extraction date:** 10/23/23 18:17:07 **Extracted by:** 152
Analysis Method : SOP.T.30.500, SOP.T.30.104.AZ, SOP.T.40.104.AZ
Analytical Batch : TE002960MYC **Reviewed On :** 10/25/23 13:44:58
Instrument Used : N/A **Batch Date :** 10/23/23 15:45:23
Analyzed Date : 10/24/23 19:12:14
Dilution : 25
Reagent : 102323.R16; 101823.R02; 101123.R02; 083123.R02; 041823.06
Consumables : 947.100; 00334958-5; 00332484-2; 1008439554; 28521042; 210823-1124; 090623; 269336; M0040957; GD220011
Pipette : TE-056 SN:21D58687; TE-060 SN:20C35457 (20-200uL); TE-108 SN:20B18337 (100-1000uL)

Aflatoxins B1, B2, G1, G2, and Ochratoxin A analysis using LC-MS/MS. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.104.AZ for sample prep, and SOP.T.40.104.AZ for analysis on ThermoScientific Atlantis TSQ with Vanquish UHPLC). Total Aflatoxins (sum of Aflatoxins B1, B2, G1, G2) must be <20µg/kg. Ochratoxin must be <20µg/kg.

Total yeast and mold testing is performed utilizing MPN and traditional culture based techniques in accordance with F.S. Rule 64ER20-39.

Metal	LOD	Units	Result	Pass / Fail	Action Level
ARSENIC	0.0030	ppm	ND	PASS	0.4
CADMIUM	0.0020	ppm	ND	PASS	0.4
MERCURY	0.0125	ppm	ND	PASS	1.2
LEAD	0.0010	ppm	ND	PASS	1

Analyzed by: 39, 30, 104, 272 **Weight:** 0.2098g **Extraction date:** 10/23/23 14:16:06 **Extracted by:** 39
Analysis Method : SOP.T.30.500, SOP.T.30.084.AZ, SOP.T.40.084.AZ
Analytical Batch : TE002953HEA **Reviewed On :** 10/24/23 14:25:17
Instrument Used : TE-051 "Metals Hood",TE-141 "Wolfgang",TE-153 "Bill",TE-157 "Bill Pump",TE-156 "Bill Chiller",TE-155 "Bill AS"
Analyzed Date : 10/24/23 11:21:47
Dilution : 50
Reagent : 050823.02; 102323.R03; 102323.R14; 102323.01; 102323.R02; 051723.06; 102323.02; 100121.01
Consumables : 12622-306CE-306C; 28521042; 210725-598-D
Pipette : TE-110 SN:20B18338 (100-1000uL); TE-169 SN: 20B16352 (Nitric Acid)

Heavy Metals screening is performed using ICP-MS (Inductively Coupled Plasma - Mass Spectrometer) which can screen down to below single digit ppb concentrations for regulated heavy metals. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.084.AZ for sample prep by microwave digestion, and SOP.T.40.084.AZ for analysis by ThermoScientific iCAP RQ ICP-MS).

	Heavy Metals	PASSED
---	---------------------	---------------

Metal	LOD	Units	Result	Pass / Fail	Action Level
ARSENIC	0.0030	ppm	ND	PASS	0.4
CADMIUM	0.0020	ppm	ND	PASS	0.4
MERCURY	0.0125	ppm	ND	PASS	1.2
LEAD	0.0010	ppm	ND	PASS	1

Analyte	LOD	Units	Result	Pass / Fail	Action Level
TOTAL AFLATOXINS	1.4870	ppb	ND	PASS	20
AFLATOXIN B1	1.4700	ppb	ND	PASS	20
AFLATOXIN B2	1.8000	ppb	ND	PASS	20
AFLATOXIN G1	1.9000	ppb	ND	PASS	20
AFLATOXIN G2	3.2500	ppb	ND	PASS	20
OCHRATOXIN A	4.6100	ppb	ND	PASS	20

Heavy Metals screening is performed using ICP-MS (Inductively Coupled Plasma - Mass Spectrometer) which can screen down to below single digit ppb concentrations for regulated heavy metals. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.084.AZ for sample prep by microwave digestion, and SOP.T.40.084.AZ for analysis by ThermoScientific iCAP RQ ICP-MS).



1231 W. Warner Road, Suite 105
 Tempe, AZ, 85284, US
 (480) 220-4470

Kaycha Labs

Social CBD Gummies Unwind Blackberry Mint 60ct 750mg

N/A

Matrix : Infused

Type: Gummy



Certificate of Analysis

PASSED

e2e Pharma

3279 E. Harbour Drive
 Phoenix, AZ, 85034, US
 Telephone: (602) 737-0077
 Email: shannon.bard@e2epharmamfg.com

Sample : TE31019007-001
 Harvest/Lot ID: 2308023B
 Batch# : 2308023B
 Sampled : 10/19/23
 Ordered : 10/19/23

Sample Size Received : 51.75 gram
 Total Amount : 1 units
 Completed : 10/25/23 Expires: 10/26/24
 Sample Method : SOP Client Method

Page 5 of 6

	Filth/Foreign Material	PASSED
--	-------------------------------	---------------

Analyte	LOD	Units	Result	P/F	Action Level
Filth and Foreign Material	0.3000 %		ND	PASS	3
Analyzed by: 93, 87, 104	Weight: 0.945g	Extraction date: 10/20/23 13:54:51		Extracted by: 87	
Analysis Method : SOP.T.40.090		Reviewed On : 10/23/23 12:15:46			
Analytical Batch : TE002936FIL		Batch Date : 10/20/23 13:38:42			
Instrument Used : N/A					
Analyzed Date : N/A					
Dilution : N/A					
Reagent : N/A					
Consumables : N/A					
Pipette : N/A					

Includes, but is not limited to: hair, insects, feces, packaging contaminants, and manufacturing waste/by-products. (Method: SOP.T.40.090 using an SH-2B/T Stereo Microscope). Testing result is for informational purposes only and cannot be used to satisfy dispensary testing requirements in R9-17-317.01(A) or labeling requirements in R9-17-317. Nor, can it be used to satisfy marijuana establishment testing requirements in R9-18-311(A) or labeling requirements in R9-18-310 - Q3.

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors.

Ariel Gonzales
 Lab Director
 State License #
 00000024LCMD66604568
 ISO 17025 Accreditation # 97164

Signature
 10/25/23

Revision: #1 - New photo



1231 W. Warner Road, Suite 105
Tempe, AZ, 85284, US
(480) 220-4470

Kaycha Labs

Social CBD Gummies Unwind Blackberry Mint 60ct 750mg

N/A

Matrix : Infused

Type: Gummy



Certificate of Analysis

PASSED

e2e Pharma

3279 E. Harbour Drive
Phoenix, AZ, 85034, US
Telephone: (602) 737-0077
Email: shannon.bard@e2epharmamfg.com

Sample : TE31019007-001

Harvest/Lot ID: 2308023B

Batch# : 2308023B

Sampled : 10/19/23

Ordered : 10/19/23

Sample Size Received : 51.75 gram

Total Amount : 1 units

Completed : 10/25/23 Expires: 10/26/24

Sample Method : SOP Client Method

Page 6 of 6

COMMENTS

* Residual TE31019007-001SOL

1 - M2 - butanes, pentanes, acetone, 2-propanol, acetonitrile, dichloromethane, hexanes, chloroform, toluene, xylenes

* Total Yeast and Mold TE31019007-001TYM

1 - Q3

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors.

Ariel Gonzales
Lab Director

State License #
00000024LCMD66604568
ISO 17025 Accreditation # 97164

Signature
10/25/23

Revision: #1 - New photo

Revision: #1 This revision supersedes any and all previous versions of this document.